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a

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
# **Certificate of Origin**

## **EC declaration of compliance**

In accordance with Appendix 11 sub A of Directive 98/37/EG which dates from 1998, we hereby declare that the listed products conforms with the published specification and complies with the requirements of the following European Directives and Standards and to their current amendments:

- Machines guideline (directive 98/37/EG which dates from 1998 as per last amendment)
- Low voltage guideline 73/23/EEG
- European standard EN 60204 and ISO 11607
- EMC guideline (directive 89/336/EEG)

## **Product**

Type

F108                  from 01-02-2000

## **Certificate of origin**

Herewith it is certified that the listed equipment has been manufactured in a member state (The Netherlands) of the European Community.

Famos bv

Paul J.M.J. van Gardingen  
General Director

## **1) In general**

All Famos rotary heat sealers are conform to the Dutch and German Directives (R3280 en DIN 58.953 Part 7) that are drawn up for sealing packing material as indicated below. All Famos F108 heat sealers can therefore be used for packing medical equipment.

The Famos heat sealers are **NOT** medical products.

It is not allowed to use the sealers in the direct presence of patients (for example, in operation theatres). The minimum distance between the heat sealing machines and patients should be 1,5m.

## **Working area**

The Famos F108 rotary heat sealer is suitable for sealing following materials:

- Laminate packing material
- Sterilisation paper packingmaterial fitted with top seal
- Tyvek (if sealer is designed as a Tyvek sealer)

For the correct temperature setting of the different packing materials we refer to the manufacturer of the materials.

The Famos F108 rotary heat sealer is **not** suitable for sealing the following material:

- Polyethylen foil
- Soft PVC foil
- Hard PVC foil
- Polyamide foil
- Polypropylen foil

(For above mentioned materials Famos bv can offer an impulse sealer).

## **Installation**

In order to guarantee a correct functioning, we advise to install the sealer at a spot:

- where the machine has enough space for ventilation
- where it is not exposed to direct sunlight or some other heatsource.
- where it is not subject to vibrations or shocks.
- that is dustfree.

Check if the on/off switch is off.

### **REMARK**


Connect the mains lead to the sealer first and than to an earthed socket with the correct Mains voltage (according to type label at the rear side of the sealer).

Switch on the on/off switch.


During heating-up the machine will mention in the display that the temperature is too low. From the off-situation onwards it takes  $\pm$  5 minutes before the pre-set temperature (at 180°) is reached.

Two bleeps indicate that the pre-set temperature is reached.  
The display will show the stand-by screen.

Feeding packages is possible now.  
The transport can be started by pressing the motor on/off switch.

#### **4) Operating panel**

##### **Explanation of the keys**

- A      Adjusting the temperature
- B      Switching on/off the transport.

##### **Operation**

###### **Display of the set temperature**

The set temperature can be seen by pressing the keys 'A' shortly.  
The set temperature will remain in the display for some seconds. After this the standard screen will appear again.

###### **Adjusting the temperature**

The temperature can be adjusted by keeping the keys 'A' pressed until the display shows the  $\blacktriangleleft$  symbol right below.  
The required temperature can be adjusted by means of the keys 'A'. When the keys are not pressed for more than 2 seconds, the adjustment will be fixed and the standard screen appears again.

### **3) Preventative maintenance**

In order to guarantee an optimal functioning of the rotary sealer and in order to prevent corrective maintenance, we advise to have the machine checked and cleaned at least twice a year. During this inspection the following parts are checked on wear and when necessary replaced.

- 12.210.300 Transportbelts
- 30.200.039 Teflon strip on heating elements
- 30.200.041 Teflon strip on guiding bars

For the article numbers and descriptions of all other parts we refer to your dealer and/or to the technical documentation.

### **Technical details**

- Temperature range	: 100 - 199 °C
- Feeder speed	: 8 m/min
- Seal width	: 15 mm.
- Self adjusting pressure	: 10 Kg (Intern mechanically adjustable: 1 – 12 Kg)
- Fuses	: See type label at the rear side of the sealer
- Capacity	: max. 365 W, average 165 W
- Main voltage	: See type label at the rear side of the sealer
- Dimensions	: 420x290x140 mm
- Weight	: 12 Kg

#### **Warning**

To avoid the risk of electric shock do not open the casing. Repairs should be carried out by a skilled and trained specialist.